# N-400

*Special microscopes*  
*Microscopi speciali*  
*Microscopios especiales*  
*Microscopes spéciaux*





# N-400FL

Epi-fluorescence microscope

Microscopi per epi-fluorescenza

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## **General features**

### The instrument

N-400FL is an epifluorescence microscope that pushes simplicity and cost-effectiveness to levels up to now considered unattainable. The basic structure is dedicated to the most demanding applications of routine fluorescence analysis. However, N400FL offers, in the same unit, brightfield and, optionally, darkfield and phase contrast capabilities, thus extending its potentials to most multi-contrast applications.

### Stand

The body is made entirely of metal, with a modern and functional design. It provides sturdy, firm and stable support under all working conditions. The microscope fits a 5-position revolver nosepiece and a large size specimen table. The focussing is performed via two coaxial macro and graduated micrometric knobs (0.002 mm / graduation) located on both sides of the stand.

### Optical system

N-400 FL is equipped with a diascope (transmitted) Kohler illuminator, and with an epi- (reflected) fluorescence excitation source. The epifluorescence optical system is implemented via the standard excitation filter-dichroic mirror-barrier filter combination, applied to a 100W Hg lamp. It is supplied with WF10x/18 mm widefield eyepieces, four plan-achromatic objectives (4x / 0.1, 10x / 0.25, 40x / 0.65 and 100x / 1.25 oil immersion), a double filterblock set (blue and green excitation). An optional a filterblock set for violet and ultraviolet excitation, and condensers and objectives for darkfield and phase contrast are available. The extensive range of optional accessories allows a quick interchange of contrast mechanisms, and it is optimized for multi-contrast observation without removal of the specimen from the stage.

### Transmitted light illuminator

The brightfield and darkfield illumination is provided by a high luminosity 6V/20W halogen bulb with intensity control. The system follows a fully-centrable Kohler scheme, for crisp and bright images. The microscope is fitted with a brightfield Abbe condenser with a numerical aperture of 1.25. Under (optional) darkfield illumination it can carry a NA1.36 condenser, thus allowing, using an optional dedicated objective, darkfield observation at 100X

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## Technical specifications

Part	Description
<b>Optical system</b>	Mechanical tube length: DIN standard - 160 mm; parfocal distance 45 mm.
<b>Head</b>	Trinocular, 30° inclined, 360° rotating. Diopter adjustment; interpupillary distance adjustment 55-75 mm.
<b>Eyepieces</b>	Wide Field WF10x/18 mm.
<b>Nosepiece</b>	5-position reversed revolving nosepiece. Ball bearing linear guides.
<b>Objectives</b>	FLUOR Plan-achromatic 4x/0.1, 10x/0.25, 40x/0.65 (spring loaded) and 100x/1.25 (spring loaded, oil imm.).
<b>Specimen stage</b>	140x132 mm 2-layer mechanical stage. X-Y travel motion: 75x55 mm with two graduated scales (0.1 mm division). Sample clip for two specimens. Coaxial controls on the right side.
<b>Focusing system</b>	Rack and pinion mechanism, with coaxial coarse and fine control knobs. Fine adjustment graduation 0.002 mm. Vertical movement range: 20 mm. Tension control on right side; upper stage drive stop on left side.
<b>Condenser</b>	Abbe condenser with double lens. N.A. 1.25. Fitted with iris diaphragm, blue filter and filter holder. Height adjustment by a rack and pinion mechanism.
<b>Illumination</b>	6V/20W with halogen bulb with brightness control. HBO 100W high pressure mercury bulb for epi-fluorescence.

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


# N-400LD




**LED FLUORESCENCE  
FLUORESCENZA LED  
FLUORESCENCIA LED  
FLUORESCENCE LED**

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## **General features**

### A revolution

Imagine a fluorescence microscope that needs a lamp change every 50 000 hours.

Imagine a fluorescence microscope with a cold light source that barely heats up during use.

Imagine a fluorescence microscope that can be switched on, used immediately, switched off and then back on again.

Imagine a fluorescence microscope that is so safe as to need no protection shield whatsoever, and that can be used by anyone, without any specific precaution.

Imagine a fluorescence microscope that can be powered by batteries, as easily as a torchlight.

Imagine a fluorescence microscope that is so sturdy and so compact that it can be used on the field, with no transport problems.

You may think that such an instrument exists in your imagination only. Actually, such microscope is real, and its name is OPTIKA N-400LD.

Developed by the OPTIKA Research labs, N-400LD marks a revolution in the field of fluorescence microscopy.

Strictly derived from model N-400FL, from which it shares the body, the optics and the filter sets, N-400LD employs high-power LED's instead of the classical mercury vapour lamp. The LED's are tailored to the specific applications (FITC-TRITC, DAPI, etc.)

The brightfield illuminator uses high-power white LED's too, and the colour temperature closely matches sunlight.

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## Technical specifications

Part	Description
<b>Optical system</b>	Mechanical tube length: DIN standard - 160 mm; parfocal distance 45 mm.
<b>Head</b>	Trinocular, 30° inclined, 360° rotating. Diopter adjustment; interpupillary distance adjustment 55-75 mm.
<b>Eyepieces</b>	Wide Field WF10x/18 mm.
<b>Nosepiece</b>	5-position reversed revolving nosepiece. Ball bearing linear guides.
<b>Objectives</b>	FLUOR Plan-achromatic 4x/0.1, 10x/0.25, 40x/0.65 (spring loaded) and 100x/1.25 (spring loaded, oil imm.).
<b>Specimen stage</b>	140x132 mm 2-layer mechanical stage. X-Y travel motion: 75x55 mm with two graduated scales (0.1 mm division). Sample clip for two specimens. Coaxial controls on the right side.
<b>Focusing system</b>	Rack and pinion mechanism, with coaxial coarse and fine control knobs. Fine adjustment graduation 0.002 mm. Vertical movement range: 20 mm. Tension control on right side; upper stage drive stop on left side.
<b>Condenser</b>	Abbe condenser with double lens. N.A. 1.25. Fitted with iris diaphragm, blue filter and filter holder. Height adjustment by a rack and pinion mechanism.
<b>Illumination</b>	High power white LED unit for transmitted light. High power LED unit for epi-fluorescence (for standard use with B and G filters).

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





**Microscopio polarizador**

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## **Technical specifications**

Part	Description
<b>Optical system</b>	Mechanical tube length: DIN standard - 160 mm; parfocal distance 45 mm.
<b>Head</b>	Trinocular, 30° inclined, 360° rotating. Diopter adjustment; interpupillary distance adjustment 55-75 mm.
<b>Eyepieces</b>	Wide Field WF10x/18 mm.
<b>Analiser</b>	Rotating filter with graduation 0°-90°.
<b>Bertrand lens</b>	Sliding in/out type with centering mechanism.
<b>Compensators</b>	1 slip (first class red), 1/4 slip, quartz wedge
<b>Nosepiece</b>	5-position reversed revolving nosepiece. Ball bearing linear guides.
<b>Objectives</b>	Strainfree achromatic 4x/0.1, 10x/0.25, 40x/0.65 (spring loaded) and 60x/0.85 (spring loaded).
<b>Specimen stage</b>	Diameter 145 mm, 360° rotating, centrable and graduated (1° division). Minimum resolution 6' when using vernier scale.
<b>Focusing system</b>	Rack and pinion mechanism, with coaxial coarse and fine control knobs. Fine adjustment graduation 0.002 mm. Vertical movement range: 20 mm. Tension control on right side; upper stage drive stop on left side.
<b>Condenser</b>	Abbe condenser with double lens. N.A. 1.25. Fitted with iris diaphragm, blue filter and filter holder. Height adjustment by a rack and pinion mechanism.
<b>Polariser</b>	Sliding in/out filter, located on the top of the illuminator.
<b>Illumination</b>	6V/20W with halogen bulb with brightness control.

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## Optional accessories

M-001	H5x eyepiece
M-020	WF10x/18mm eyepiece
M-006	EWF10x/20mm eyepiece
M-003	WF16x/12mm eyepiece
M-004	Micrometer eyepiece WF10x/18mm
M-005	26x76 mm micrometric slide. Range 1 mm, div. 0,01 mm
M-410	4x fluorescence objective
M-411	10x fluorescence objective
M-605	20x plan achromatic objective
M-412	40x fluorescence objective
M-689	60x plan achromatic objective
M-413	100x fluorescence objective
M-050	Polarising set (filters only)
M-051	Rotating table for polarising set
M-057	Darkfield condenser for dry objectives
M-058	Darkfield condenser for immersion objectives
M-059	100x plan objective with iris diaphragm for darkfield purpose
M-045	Simple phase contrast set 10x
M-046	Simple phase contrast set 40x
M-047	Simple phase contrast set 100x
M-048	Complete phase contrast set with plan obj. 10x, 20x, 40x, 100x
M-062	Photo tube adapter for SRL cameras
M-063	CCD camera adapter
M-150	Additional filter set U (ultra-violet) and V (violet)
M-153	V (violet) LED illuminator 385 nm for N-400LD
M-154	UV (ultra-violet) illuminator 365 nm for N-400LD
M-014	6V/20W halogen bulb
M-151	HBO100W high-pressure mercury bulb for fluorescence
M-032	Dust cover type 4
M-974	Blu filter, 32 mm diameter
M-976	Green filter, 32 mm diameter
M-978	Yellow filter, 32 mm diameter
M-988	Frosted glass filter, 32 mm diameter

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